

The logo for 'evo' is displayed in white lowercase letters on an orange curved background.

# Inogen One G2

Oxygen System (E1390 & E1392)

IS-200

First there was stationary. Then there was portable. Now there is the Inogen One G2, the Independent System that defines a new category of oxygen device. With the ability to perform as both a stationary and a portable device, the G2 is engineered to meet the changing demands of the oxygen industry.



FAA  
Approved

- **Small** — 42% smaller than the Inogen One.
- **Lightweight** — only 7 lbs. with the battery.
- **Performance** — delivers more oxygen 900 mls, faster and at the most advantageous time for absorption.
- **Reliability** — built to last 25,000 compressor hours.
- **Battery** — 4 hour or 8\* hour battery options.
- **Virtually Silent** — operates so quietly, you hardly know it is running.



## Inogen One G2 Specifications

Model:	IS-200 Inogen One G2 System
FDA#:	K032818
Weight:	7.0 lbs with battery
Flow:	900 ml/min
Warranty:	3-Years on device 1-Year on battery
Size:	10.1cm L x 27.4cm W x 24.1cm H (40 mm x 105 mm x 93 mm)
Battery:	12 cell battery provides 4 hours 24 cell battery provides 8 hours*



## Inogen One G2 Includes

The Inogen One G2 System includes the concentrator, carry bag, universal power supply, a 12 cell battery, cart, nasal cannula and user manual.

## Why the Inogen One G2?

### **For Patients—**

the Inogen One G2 allows for a **more convenient** and **more active life** both at home and on the road.

### **For Providers—**

the Inogen One G2 means **improved business efficiency.**

## Inogen One vs Inogen One G2

Size (cm)	29.5 x 15.2 x 31.5	27.4 x 10.1 x 24.1
Volume	9000 cc	5200 cc
Weight	10 lbs	7.0 lbs
Noise	<40 dba	<38 dba
12 cell battery	3 hr / 2 hr	4 hr / 2 hr
24 cell battery*	NA	8 hr / 4 hr
Oxygen	750 ml/min	900 ml/min.
Conserver	Night enabled	Night enabled
Data Log	2-3 months / Infrared	1 year / USB or RF
Device Warranty		3 year
Battery Warranty		1 year



**evo** Medical Solutions

2636 289th Place  
Adel, IA 50003  
515.993.5001  
515.993.4172 f  
800.759.3038 (toll free)  
www.evomedical.com

\*24 cell battery available Spring 2010.

InogenOne vs. Inogen One G2 device specifications and information provided from Inogen.

P/N 22435 Rev A